



SPECIFICATIONS

- Handle:** Plastic, with media indication according to EN13792
- Spout:** Full swivel action
- Hose nozzle:** Fixed plastic
- Headwork:** Rubber (wrist action lever: ceramic)
- Opening/closing function:** 2 x 360° (wrist action lever: 90°)
- Temperature range:** 0-65°C
- Max. working pressure:** 10 bar / 147 psi
- Test pressure:** 1 x working pressure
- Weight:** 1.4 kg

MATERIALS WITH MEDIA CONTACT

- Metals:** Brass
- Plastics:** PP, PTFE, SPX (hose)
- Sealing:** EPDM
- Lubricant:** Silicone based
- Ceramics:** Al₂O₃ 96% (wrist action lever only)

ORDERING INFORMATION

For standard configuration please refer to the item numbers below.

Other options:

- aerator
- removable plastic hose nozzle
- integrated non-return valve
- fixed spout (rear, left or right)
- alternative hose inlet connection / without a hose
- other inlet options (type, length)
- national approvals
- WPH (Water Potable Hot)
- WNH (Water Non-Potable Hot)

BROEN UNIFLEX™ FITTINGS FOR WATER

With wrist action lever



Model: Plastic handle



Media	S = 150mm	S = 200mm	S = 250mm	S = 300mm
Potable water, cold	2534 008 1001-52B	2534 084 1001-52B	2534 085 1001-52B	2534 086 1001-52B
Non-potable water, cold	2534 008 1007-52B	2534 084 1007-52B	2534 085 1007-52B	2534 086 1007-52B

Model: Wrist action lever



Media	S = 150mm	S = 200mm	S = 250mm	S = 300mm
Potable water, cold	2534 087 1001-52B	2534 088 1001-52B	2534 089 1001-52B	2534 083 1001-52B
Non-potable water, cold	2534 087 1007-52B	2534 088 1007-52B	2534 089 1007-52B	2534 083 1007-52B

For further information and mounting instructions please refer to www.BROEN.com or contact your local BROEN LAB representative.