

Laboratory PU Lab Chair ILC01



Product Description

- ~Colour : Black / Light Grey
- ~Back : PU back
- ~Seat : PU seat
- ~Footrest : Nylon footring in diameter 450mm
- ~Mechanism : One lever sincro 2.0mm thickness mechanism ith chromed bracket
- ~Gaslift : 260mm chromed gaslift
- ~Base : 300mm chromed base
- ~Caster : 50mm nylon caster
- ~Net weight : 10.7 kg/ctn
- ~Gross weight : 11.7 kg/ctn

Product Dimension

W465mm X D530mm X H970mm~1210mm

Back size : 420mm X 310mm

Seat size : 470mm X 450mm

