

DILUTION TANK



	Tech. Drawing	Dimension		
		Cat. No.	PP110-010	PP110-020
		A	230mm	230mm
		B	121mm	121mm
		C	318mm	318mm
		D	244mm	244mm
		E	168mm	168mm
		F	51mm x 1	51mm x 3
		G	38mm x 2	—

		Dim.	38mm		38mm Flexible Connector 0.3m 0.4m 0.6m 1.0m (It is ideal for underbench and mobile fumecupboard use.)	PP080-004 PP080-003 PP080-001 PP080-002
		A	0.3m	0.4m		
			0.6m	1.0m		

		Dim.	38mm	51mm	FI to Pipe Coupler 38mm 51mm	PP325-010 PP325-020
		A	1 1/2" BSP	2" BSP		
		B	6	6		