

HAP FUME HOOD

- Work Surface: 10 mm high anti-corrosive PP panel.
- Main Structure: 8 mm high anti-corrosive PP panel.
- Window Glass: 6 mm transparent tempered glass.
- Exhaust Canopy: hydrodynamic design canopy with PP material.
- Cabinet Door: 12 mm high anti-corrosive PP panel with PVC handle.
- High anti-corrosive plastic hinge.
- Water in and drain components:
 - a. Cock: brass alloy with powder coating.
 - b. Control Valve: brass alloy with powder coating.
 - c. Sink Cup: PP material.
 - d. Drain Trap: PP material.
- Control Panel:
 - a. Blower and lamp switches
 - b. 230V electrical outlets.
 - c. Adjustable Overload protection and magnetic switch set.
 - d. Blower and electrical outlet breakers



| Model | Dimensions |
|---------|---------------------|
| HAP-090 | W900*D800*H2350 mm |
| HAP-120 | W1200*D800*H2350 mm |
| HAP-150 | W1500*D800*H2350 mm |
| HAP-180 | W1800*D800*H2350 mm |