

---

## 944ml Saline Eyewash Bottles with Wall-Mount & Sign

SKU: 944ml Saline Eyewash Bottles with Wall-Mount & Sign (EWB944ECS)



### Product Information

---

## Categories: Personal Eye/Facewash

EWB944ECS personal emergency eyewash bottles are designed as a secondary eye wash device and is intended for support use, not as a replacement for self contained or plumbed units.

### Product Description

---

[vc\_row][vc\_column width="1/2"][vc\_column\_text]The EWB922ECS saline eyewash squeeze bottle should be located in the immediate vicinity of employees working in a potentially hazardous area. It is used to deliver initial first aid while getting the person to a primary unit, or as extended aid when transporting a person from the primary eye wash unit to a medical facility.

Time is critical, the range of eye wash bottles have been designed to ensure that opening and commencement of flushing is achieved in optimal time. Just twist the top to break the seal and commence flushing.

Contents:

1 x 944ml saline soft squeeze bottle

1 x wall mount

1 x eyewash sign[/vc\_column\_text][vc\_column][vc\_column width="1/2"][/vc\_column][vc\_row]