

Table Top Low Speed Centrifuge



Technical Parameters:

Model	BKC-TL5E
Max.Speed	5000rpm
Max.RCF	3500×g
Max.Capacity	4*2*100ml
Speed Precision	±1% or ±20rpm
Time Range	1s~999min
Power Supply	Standard: AC220V 50Hz Optional: AC110V 60Hz
Consumption	500W
Noise	≤65dB(A)
External Size(L*W*H)	540*410*330mm
Package Size(L*W*H)	700*530*480mm
Net/Gross Weight	40kg/55kg

Rotor:

Type	NO	Capacity	Max.Speed(rpm)	Max.RCF(×g)
Swing Rotor	NO.1	4*1*50ml	4000	2846
	NO.2	4*1*100ml		3007
	NO.3	4*2*50ml		2846
	NO.4	4*2*100ml		3007
	NO.5	4*4*2/5/10/15ml		2705
	NO.6	4*6*2/5/10/15ml		2705
	NO.7	4*8*2/5/10ml		2705
Angle Rotor	NO.8	8*15ml	5000	3300
	NO.9	12*10ml		3440
	NO.10	12*15ml		3440
	NO.11	12*20ml		3360
	NO.12	18*10ml		3500
	NO.13	24*10ml		3500
	NO.14	4*50ml		3300
	NO.15	6*50ml		3300
	NO.16	4*100ml		3300

Features:

- * Overspeed, unbalanced protection, safe and reliable.
- * Large torque motor direct drive, no carbon pollution, maintenance free, long service life.
- * OLED digital display, real-time display of all operating parameters.
- * Embedded microprocessor control, high control precision.
- * 10 levels speed-up and can be adjusted within a certain range.
- * Automatically memory the last parameters and restore the parameters before shutdown.